


MDO-FB Front-Basis Distractor

Non-extraction therapy for anterior crowding of the mandible


Developed in cooperation with Albino Triaca, MD, DMD and Roger Minoretti, MD, DMD,
Center for Maxillofacial Surgery, Pyramide Clinic, Zürich, Switzerland



- Resolution of mandibular anterior crowding without extraction – no negative impact on the facial profile
- Correction of unfavorable relations between the anterior dento-alveolar area and the skeletal base
- Combination of translation and rotation to achieve an optimal position and axis of the anterior mandibular dentition
- Application in local anesthesia

Orthognathics MDO-F – Front-Basis Distractor

The Orthognathics MDO-FB Front-Basis Distractor allows a correction of position and axis of the anterior mandibular dentition



dento-alveolar area and the skeletal base. The device allows a combination of translation and rotation to achieve an optimal position and axis orientation for the anterior mandibular dentition.

MDO-FB is designed for application in the mandible with retro-alveolism, for example skeletal class I with a dental class II, where a gap is created in the area of the premolars and where a deep sub-labial fold is present which needs to be smoothed out by the advancement of the base of the anterior alveolar ridge.

The MDO-FB Front-Basis Distractor allows a gradual repositioning of an anterior dentate segment of the mandible using distraction osteogenesis. It facilitates a non-extraction therapy for patients with anterior tooth crowding in the mandible or the correction of an unfavorable relation between the anterior

Technical Data:

Material	surgical steel
Bone screws	for 2 mm screws
Transport screw	2 mm screw

For any questions on this product, please send us an e-mail or write us by fax or mail.

- ↑ Send more information
- ↑ Send price quotation
- ↑ Please arrange a meeting or demonstration/surgery
- ↑ Please keep me informed about workshops and presentations on MDO-H and other Orthognathics products

Orthognathics Ltd.
 Schuetzenstr. 41
 CH-8304 Wallisellen/Switzerland
 Tel. +41 (0)43 2689200
 Fax +41 (0)43 2689202
 info@orthognathics.com

Name: _____
 Address: _____

 Phone: _____
 Fax: _____
 e-mail: _____

 Date Signature