

MDO-O Orthodontic Anchor

Maximal anchorage in the mandible

Developed in cooperation with Albino Triaca, MD, DMD, Roger Minoretti, MD, DMD and Michele Antonini, DMD, Center for Maxillofacial Surgery, Pyramide Clinic, Zürich, Switzerland



- Absolute anchorage for orthodontic treatment in the mandible
- Provides anchorage in difficult cases, e.g. with periodontitis or missing teeth
- Simple placement, short chair time, versatile application

Orthognathics MDO-O – Orthodontic Anchor

The Orthognathics MDO-O Orthodontic Anchor allows bone-anchored orthodontic movements in the mandible



The MDO-O Mandibular Orthodontic Anchor provides the necessary absolute anchorage to accomplish larger orthodontic tooth movements without the danger of side-effects on anchorage teeth or the requirement of extra-

oral anchorage aids. It therefore allows predictable treatment also in difficult cases, e.g. with periodontal problems or missing teeth.

The MDO-O can be placed in a small, minimal invasive procedure and the loading can be started within a few days. MDO-O is designed to work with standard metal brackets.

Technical Data:

Thickness	1 mm
Material	surgical steel
Bone screws	for 2 mm screws

For any questions on this product, please send us an e-mail or write us by fax or mail.

- ? Send more information
- ? Send price quotation
- ? Please arrange a meeting or demonstration/surgery
- ? Please keep me informed about workshops and presentations on MDO-O and other Orthognathics products

Orthognathics Ltd.
 Winterthurerstr. 92
 CH-8006 Zurich/Switzerland
 Tel. +41 (0)43 2689200
 Fax +41 (0)43 2689202
 info@orthognathics.com

Name: _____

Address: _____

Phone: _____

Fax: _____

e-mail: _____

Date

Signature