

MDO-R

Rapid Palatal Expansion Distractor

Precise Palatal Expansion


Developed in cooperation with Albino Triaca, MD, DMD, Roger Minoretti, MD, DMD and Michele Antonini, DMD, Center for Maxillofacial Surgery, Pyramide Clinic, Zürich, Switzerland



- Precise, bone-anchored rapid palatal expansion
- Minimal invasive design
- Expandable from initially small size to a long expansion distance

Orthognathics MDO-R – RPE Distractor

The Orthognathics MDO-R Rapid Palatal Expansion Distractor allows bone-anchored palatal expansion in the Maxilla



The special pin anchorage on the hard palate allows for a minimal invasive insertion and avoids any unnecessary irritation of the palatal mucosa during the treatment. The design of the device enables a placement even in very narrow palates while at the same time long distraction distances are possible without any changes in the device.

The MDO-R Rapid Palatal Expansion Distractor allows a fast and precise widening of the Maxilla. It avoids any negative side effects on the maxillary dentition as usually found with the typical dental based appliances.

Technical Data:

Distraction distance	11+ mm
Material	surgical steel

For any questions on this product, please send us an e-mail or write us by fax or mail.

- ? Send more information
- ? Send price quotation
- ? Please arrange a meeting or demonstration/surgery
- ? Please keep me informed about workshops and presentations on MDO-R and other Orthognathics products

Orthognathics Ltd.
 Winterthurerstr. 92
 CH-8006 Zurich/Switzerland
 Tel. +41 (0)43 2689200
 Fax +41 (0)43 2689202
 info@orthognathics.com

Name: _____

Address: _____

Phone: _____

Fax: _____

e-mail: _____

 Date Signature